

RBM7-RNA binding motif protein 7 Gene

Reference: PTXBC034381

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNA-binding protein 7; RNA binding motif protein 7		
For research use only.			

GENE DESCRIPTION	RNA binding motif protein 7
GENE ID	10179
ORIGIN SPECIES	Human
ACCESSIONS	BC034381
VECTOR	pDONR223
NCBI SYMBOL	RBM7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcgggcgggcggaagcggatcgactctcttgggcaacctgaaacgaaagtacgaggagctctttcagctttccaccagctggccagtaataaagggtgaaattccaaagataaggatggtaaaccaagcagttgcgtttggaattcaaacatgaagt
gtctgttcttatgcaatgaatctacttaattgaatcaaaacttatggaagccctatcaaaatcaatttagatcagggaagttagtcatgcccaagaatgtcagtttgcataatcccaacatcatgttggaattcaagccctacctccacatctctagcaggtaacgaaaggactatggataacat
gacttcacgacagataaattcagagatcttctctccagaaaatttcagagacaaagcagtgatgaacagtgctttgagacaaatgtcatatgggtggaatttgggtcttcacctctggatcaatcaggattttaccatcagttcaatcacagtcagtttcaatcagttctcaagctccc
agtggcgccaaggtacacacatcatcacagcgtaaagtcagaatgaattcttaccctacctagcagatagacattatagccgggaacagcgttacactgatcatgggtctgaccatcattacagaggaaagagagatgatttctatgaagacaggaatcatgatgactggagccatgacta
tgataacagaagagacagtagtagagatggaaatggcgctcatctcgacactaa