

FAM98C-family with sequence similarity 98, member C Gene

Reference: PTXBC034383

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM98C; family with sequence similarity 98 member C		
For research use only.			

GENE DESCRIPTION	family with sequence similarity 98, member C
GENE ID	147965
ORIGIN SPECIES	Human
ACCESSIONS	BC034383
VECTOR	pDONR223
NCBI SYMBOL	FAM98C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtccaagaactggacctaccctccaagccctgggctgccagacctgacccaggacccccgccagccagctgctgaggttgcatgctaagatctcagagctgcagcttctgccccagggtccctgcagccctctcagctgctcgtatgatcacccagatggga
agcgttggagttctgttccaaagcctcagagatcagtagctgctgcccgcctgctcctcctcaagccttgacctactacatctgcttccactggagtgaccgggcagagggcccaaggagagggccatgaggcagtgctgatcccaattcagaggttctgacccagaatcgga
catctcattgacacagttctggctgcccagccgacctgtccgagaggacctgtgtgcatcaacaagggtcttatgggcaacgttccagaccggggggccgcccgaatgagctggagcctccatgccacctggaggagccgaagagaggtggagggcccgatgttgg
ggctgcaagaagaagaagaagtaaggggactggtggtcggggcggggggtcctcatgatgctaa