

RNF175-ring finger protein 175 Gene

Reference: PTXBC034385

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RING finger protein 175		
For research use only.			

GENE DESCRIPTION	ring finger protein 175
GENE ID	285533
ORIGIN SPECIES	Human
ACCESSIONS	BC034385
VECTOR	pDONR223
NCBI SYMBOL	RNF175
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgggggacggcgccggcgaaggcagcgggtgctggaggccccccgcagcaggagcagctctctacaaaagcttctgcagaagacacatggaacctgcagcaggagagggatgtacaagatgcacggggccacgattccatgcacgtggaatgatcttgatctt
cctctgcgttctggtcattgccagatagtgctggttcagtgagacagaggcatggccgatcctacaatcgtgtgacctgttgcagatgtgggtgtgctcccttataattcagataaaattactggtggcggtttctgtctatgtgggggatgttctccgttattaccagttacatcctctcaga
gctaccggaaaacccctctcagggaaggacaccagattggtctacaatggttctttgatctacaactcagctatgcatgttggtgtgtgtggttactggcgatcatgtttacaatgtgtggttcaatctgtttttcaaaatcaaaagtagagattccatggatttggcattgtctttgttctacg
gcctctactatggagtaatggggagagactttccgagatctgctcagactacatggcttccactataggggttctacaggtgcagccgggttgctacaaggagcttatcgacaatctgtgcagctctgtgggcagaagatcattgtggagcttgatgaagaagggtcattgaaacaccta
ccagcttctgtaatcatgtcttcatgaattctgcacccgaggtgtgtgtatcgttgggaaaaagcagacttgccttactgcaagagagaagtgattgaagaggatgatcagtaacccctgggagcgacacatttctgtatggacaatcctggattggcttctgtatttggcctggca
acctgtggtgataggaatagttcaaggcattactatcactagggtggaatag