



Reference: PTXBC034403

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD182; CD183; CKR-L2; CMKAR3; GPR9; IP10-R; Mig-R; MigR; C-X-C chemokine receptor type 3; G protein-coupled receptor 9; IP-10 receptor; Mig receptor; chemokine (C-X-C motif) receptor 3; chemokine receptor 3; interferon-inducible protein 10 receptor; C-X-C		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) receptor 3
GENE ID	2833
ORIGIN SPECIES	Human
ACCESSIONS	BC034403
VECTOR	pDONR223
NCBI SYMBOL	CXCR3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]