

PDCL2-phosducin-like 2 Gene

Reference: PTXBC034431

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GCPHLP; phosducin-like protein 2; phosducin like 2		
For research use only.			

GENE DESCRIPTION	phosducin-like 2
GENE ID	132954
ORIGIN SPECIES	Human
ACCESSIONS	BC034431
VECTOR	pDONR223
NCBI SYMBOL	PDCL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggatcccaatgaagatacagaatggaatgacattttaagagatttcggcattctctcctaaagaagagtcacaaagatgaaattgaagaatgggttttacgaaagaagcaatgggtaaacatttgaaaagatgactcttcacagctaaaggagctgaagatgaattaat
gaagaagatatgcaggctgttgaacatatagaagaagcgttacagggaatggaaagctcttaagaaaaacaaaaatttgagaattaagagaaattctggaaatcagatgtgaatgaagtcacaaatgcagaagaagatgtgtgggtataattcatctacagatcaagcatcccaa
tgtgtttgtgttaaccagcatcttagtctctagcaagaagttccagaactaaattgttaagccatcgtgaatagctgtattcaacactaccatgacaattgtttaccaacaattttgtgtataaaaatggtcagatagaagccaaatcattggaattatagaatgtggaggataaatctca
agctggaagaacttgaaatggaagctagcagaagttggagcaatacagactgatttgaagaaaacccagaaaagacatggtagatatgtgtatctcaattagaacactctattcatgatgacagtgatagctccaacagtgataatgaaccaatagagagaatattcaataatag
cttttag