



Reference: PTXBC034435

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZC3HDC3; zinc finger CCCH domain-containing protein 3; smad-interacting CPSF-like factor; zinc finger CCCH-type domain containing 3; zinc finger CCCH-type containing 3		

For research use only.

GENE DESCRIPTION	zinc finger CCCH-type containing 3
GENE ID	23144
ORIGIN SPECIES	Human
ACCESSIONS	BC034435
VECTOR	pDONR223
NCBI SYMBOL	ZC3H3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]