



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SPANXB2-SPANX family, member B2 Gene

Reference: PTXBC034472

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

For research use only.

GENE DESCRIPTION	SPANX family, member B2
GENE ID	100133171
ORIGIN SPECIES	Human
ACCESSIONS	BC034472
VECTOR	pDONR223
NCBI SYMBOL	SPANXB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccaacaatcagctgctccaggctgaagaggagctgccctgtgaatccaacgaggccaacgaggccaatgaggccaacaagcagatgccggagaccccaactggggactcagaccgcaactctcctaaaaaagaaaaacatctgctctgaccatactagtgttcctacagagggaacgtgaaaaaacaatctccagaagaaactgtgtaatgaccacgccgagagaacaagaatcaacccccaccaaatgagagaggaagaaatcatagaataaacgactgaagaagacacaaaaagtag