



Reference: PTXBC034488

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABC27; ATP-binding cassette sub-family F member 1; ATP-binding cassette 50 (TNF-alpha stimulated); ATP-binding cassette, sub-family F (GCN20), member 1; TNF-alpha-stimulated ABC protein; TNFalpha-inducible ATP-binding protein; ATP binding cassette subfami		

For research use only.

GENE DESCRIPTION	ATP-binding cassette, sub-family F (GCN20), member 1
GENE ID	23
ORIGIN SPECIES	Human
ACCESSIONS	BC034488
VECTOR	pDONR223
NCBI SYMBOL	ABCF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]