

FAM124A-family with sequence similarity 124A Gene

Reference: PTXBC034497

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM124A; family with sequence similarity 124A; family with sequence similarity 124 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 124A
GENE ID	220108
ORIGIN SPECIES	Human
ACCESSIONS	BC034497
VECTOR	pDONR223
NCBI SYMBOL	FAM124A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccacaaagcggggcggcgaggagacgactgcgtggactggcgccgagacggagggtccgactacagccactgtcctccagcagcagtgagcttccgttgaagagcgcaggacccttctggtcagcatccacataatcgagaccaggaggagtc
agccccctgaggaggccatcgacaacgtcctggcgtggatccaccccgacccctggttccgggtgctcgagagggcggtcccgcgggcggaagcccccaaggcgctcagccagcgtggtggtgctgttctgcaggaggagtacggcgaagcagatcct
gcagctgcaccgcacactgcagcagccgccccggcgccaccacacacggagcaggtgcacggcggttctgccccctacgtccctgcagccaggacttctacgctggccccgggagccgcttgggcatccggccccgtgactacggcgaagaaatcgtgccttaccgt
ctactgtcgtacgacaactatgctgacagcctcaggttctaccagctgattctccggaggagccccagccagaagaagcggacttctgcatcttccctatttttccaacctggatgtggacatccagttctccctgaaaagactgcctgtgaccagtgcctgggtgccaccgactcctcc
gtgctggagttccgagtgaggacatagggcagctgctcctcctgccccaccttgcagccccatcagcagggcgctggcagacggaggaccatgatgggaacaagatcctctacaggtactgggggagccgtgtctgtctgtttag