

ABHD12B-abhydrolase domain containing 12B Gene

Reference: PTXBC034603

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein ABHD12B; BEM46L3; C14orf29; c14_5314; abhydrolase domain-containing protein 12B; alpha/beta hydrolase domain-containing protein 12B; abhydrolase domain containing 12B		

For research use only.

GENE DESCRIPTION	abhydrolase domain containing 12B
GENE ID	145447
ORIGIN SPECIES	Human
ACCESSIONS	BC034603
VECTOR	pDONR223
NCBI SYMBOL	ABHD12B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgggtgcaagtatcaattatccctgttaaaagtaccggaacattccaggattttacgtacattatggatgccctgagaaaagacaaaataatcttctaatgatgaaatgttaaattccctttctctctctctcatctacatggaggatgacaggacagtccttggagtatgg
aaaaagctctatgaaattgcacgaatgcatacaggaacaaagagagggtcaagatggttatcttctctctggttccaacacacacactgcttftaaaagccccacactgttaataaccgtgagagatttcttgagcaagcagtggtcatga