

S100A7-S100 calcium binding protein A7 Gene

Reference: PTXBC034687

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PSOR1; S100A7c; protein S100-A7; psoriasin 1; S100 calcium binding protein A7		
For research use only.			

GENE DESCRIPTION	S100 calcium binding protein A7
GENE ID	6278
ORIGIN SPECIES	Human
ACCESSIONS	BC034687
VECTOR	pDONR223
NCBI SYMBOL	S100A7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcaacactcaagctgagaggtccataatagcatgatcgacatgtttcacaatacaccagacgtgatgacaagattgacaagccaagcctgctgacgatgatgaaggagaacttcccaactccttagtgcctgtgacaaaaaggccacaattacctcgccgatgtttgagaa
aaaggacaagaatgaggataagaagattgattttctgagtttctgtccttgcctggagacatagccacagactaccacaagcagagccatggagcagcgccctgttccggggcagccagtga