

ULBP2-UL16 binding protein 2 Gene

Reference: PTXBC034689

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALCAN-alpha; N2DL2; NKG2DL2; RAET1H; NKG2D ligand 2; retinoic acid early transcript 1H; UL16 binding protein 2		

For research use only.

GENE DESCRIPTION	UL16 binding protein 2
GENE ID	80328
ORIGIN SPECIES	Human
ACCESSIONS	BC034689
VECTOR	pDONR223
NCBI SYMBOL	ULBP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcagccgctaccaagatcctctgtgctcccgcttctgtctcctgtctgctccgctgtgctccgctggcgagccgacctcactctctttgtatgacatcacccgtcatccctaagtcagacctggaccaggtgtgtgctggtcaagccaggtggatgaaaagactttc
ttcactatgactgtggcaacaagacagtcacacctgtcagtcctcctggggaagaactaaatgtcacaacggcctggaaagcacagaaccagctactgagagaggtgtgtgacatactacagagcaactgcgtgacatcagctggagaattacacccaaggaaacctcaccctgc
aggccaggtgtctgtgagcagaagctgaaggacacagcagtgatcttggcagttcagtttcgatggcgagatcttctcctcttgaactcagagaagagaatgtggacaacggttcctctggagccagaagatgaaagaaaagtgaggagaatgacaaggtgtggccatgtccttc
cattacttcaatgggagactgtataggatggcttgaggacttcttgatggcagtgacagcaccctggagccaagtgcaggagcaccactgcatgtctcagcacaaccaactcaggccacagccaccctcatcttctgctcctcatctcctcctgttcctcc
ctggcatctga