



Reference: PTXBC034703

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EDG8; Edg-8; S1P5; SPPR-1; SPPR-2; sphingosine 1-phosphate receptor 5; S1P receptor 5; S1P receptor Edg-8; endothelial differentiation G-protein-coupled receptor 8; endothelial differentiation, sphingolipid G-protein-coupled receptor, 8; sphingosine 1-ph		

For research use only.

GENE DESCRIPTION	sphingosine-1-phosphate receptor 5
GENE ID	53637
ORIGIN SPECIES	Human
ACCESSIONS	BC034703
VECTOR	pDONR223
NCBI SYMBOL	S1PR5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr