



Reference: PTXBC034709

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CX30.3; EKV; gap junction beta-4 protein; connexin 30.3; gap junction protein, beta 4, 30.3kDa; gap junction protein beta 4		

For research use only.

GENE DESCRIPTION	gap junction protein, beta 4, 30.3kDa
GENE ID	127534
ORIGIN SPECIES	Human
ACCESSIONS	BC034709
VECTOR	pDONR223
NCBI SYMBOL	GJB4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]