



Reference: PTXBC034761

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	7TMR; ADMR; AM-R; AMR; G10D; gamrh; hrhAMR; G-protein coupled receptor 182; adrenomedullin receptor; G protein-coupled receptor 182		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 182
GENE ID	11318
ORIGIN SPECIES	Human
ACCESSIONS	BC034761
VECTOR	pDONR223
NCBI SYMBOL	GPR182
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1