

TRPV6-transient receptor potential cation channel, subfamily V, member 6 Gene

Reference: PTXBC034814

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABP/ZF; CAT1; CATL; ECAC2; HSA277909; LP6728; ZFAB; transient receptor potential cation channel subfamily V member 6; Alu-binding protein with zinc finger domain; calcium transport protein 1; calcium transporter-like protein; epithelial apical membrane ca		

For research use only.

GENE DESCRIPTION	transient receptor potential cation channel, subfamily V, member 6
GENE ID	55503
ORIGIN SPECIES	Human
ACCESSIONS	BC034814
VECTOR	pDONR223
NCBI SYMBOL	TRPV6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtacgttgagaatcactgtccaggcctgcattactcctcagctctggggcagagggaagccagcccaagcacgggctggcaggcgtgaggaactctcctgtggcctgtcatcaccttccgacaggagcactgcatgtcagagcacttataaacaggccagcctgtgggcgtcgggtctccacccagggtcataagtgaggagagagccctccagggcaccagcagggtgcagggaagtgcagagcttgggaagcgtgtgagtgaggagacagggaacgctctgggggtgggaagtggggctag