

MGC34796-SPR pseudogene Gene

Reference: PTXBC034822

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	SPR pseudogene
GENE ID	414927
ORIGIN SPECIES	Human
ACCESSIONS	BC034822
VECTOR	pDONR223
NCBI SYMBOL	MGC34796
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgctgcggcagctggagccgagctgggcccggcggtgacttgcgcgtgtgcgggtgcccgagacctggcgccgagctggcttcagcccgctgttggggccttgcgtcagctcccaggaccaggagctgcagcaactgctctcatcaacaacgcccgtc
tcctggggatgtgtccaaaggcttcgtggacctgggtgactccactcaagtgaacaagactggcgctgaactgacctccatgctctcttgcgggacttcagcgtctcgaaggccttaccaccagctcctggcctcaacagaactgtgttaa