

## CCPG1-cell cycle progression 1 Gene

Reference: PTXBC034914

|                |                                                                                                                |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                            | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CPR8; cell cycle progression protein 1; cell cycle progression restoration protein 8; cell cycle progression 1 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | cell cycle progression 1                                                                         |
| <b>GENE ID</b>          | 9236                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC034914                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | CCPG1                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgtcaacttttttctaagtattttctcagcctagtgacgatgatcaagtagtgatgaacagtagtcaagcagctccttttagacgacgccgtgctaggaagaagaccgtttctgctcagaatctgaagaccggctagttgctgaacaagaactgaaccttctaaaggattgag  
taaacgtcagttcagtagtgggtctcaataagtggttatacttgccttggtgattgcaatcagcatgggatttggcatttctatggcacaattcagattcagaagcgtcaacagttagtcagaaagatacatgaagatgaattgaatgatataagggattatcttccagtgtaacagggaacaaga  
atcttttatagattataagtcattgaaagaaatcttgcagggtgttgacacttactgaagcagagagaagatgtcctttgaaactcagaaacgaaccttgcacagaaatcagtttaagagtatccctggagaagggaagaaagccctatctcattacaggaagagttaaacaaactaag  
agaacagattagaataattggagataaagggaacaaagctgaattagtttaagaaatcagaaacttaagcagcatttgggaagggaaagcagaaaaaacacagctttctagtcagaaggagactctgttgacagaagcaaatgatctaaagagagaactggagagagaacgactagt  
aactacggctttaagggggaactccagcagtttaagtggtagtcagttacatggcagtcagatttcccaatgtatatactgaaaaaaaggaaatagcaactcttacgggaagactcactgagctggaactgaagtaacctcgaacagcagcgttctgatttgggaaagattgtatgtt  
aggcaaaagatcaaatggaaaacaggaacagatggaaaaaagaaggggcagagggaagccacagggcctaaaaataagtcagaaagaaacatttttgggtcagtttaaggaaacatttgatgccatgaagaattctaccaaggagtttgaaggcatataagagaaaattaaagcag  
gctaaagaagctgtgaaggaaaatctgaaaaatctcagattcagtttaaatccacttcagacactttaagataccaccagaataatcttggatgaaaagggttaataaaagatttgggtctacaaaagaagcagctgaaaaaccaagaacagtttttagtgactatttaccacagataaagg  
cacctacagaaacattataatagagccctactatgcaaaatgatggaaggaagaaagccagttcactttaagaattcagaaaaatacaaatcaaaagaatgcagtcctggcagattgtagagaaaattctcattcttcagagaggcgttctgtgtgttattgtgtcaca  
agagtcagagccttttaacatagtggtgaatcctataaggatggatgaatttagacagataattcaagggtacattttaaagaactggatacttttgcactgggaacgaactgacagttcatcaataagttttcctaaacgggtgtctttatcatgacagaagctctcactgacttgtta  
gatgttaagaattatcttagaacaatgaaggaatatgaagtagataatgatggagatttggagaagttggatgaatatatatagacacttcttggtcacactttttccctccatagaccaggtcgggttacataaaacccgtgtcattacagtagttgtaa