



Reference: PTXBC034948

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NT5; orphan sodium- and chloride-dependent neurotransmitter transporter NTT5; CTC-301O7.4; solute carrier family 6 (neurotransmitter transporter), member 16; solute carrier family 6 member 16		

For research use only.

GENE DESCRIPTION	solute carrier family 6, member 16
GENE ID	28968
ORIGIN SPECIES	Human
ACCESSIONS	BC034948
VECTOR	pDONR223
NCBI SYMBOL	SLC6A16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]