

C1orf49-chromosome 1 open reading frame 49 Gene

Reference: PTXBC034972

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf49; TSC24; testis-expressed sequence 35 protein; Testis-Specific Conserved gene 24kDa; testis expressed 35		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 49
GENE ID	84066
ORIGIN SPECIES	Human
ACCESSIONS	BC034972
VECTOR	pDONR223
NCBI SYMBOL	C1orf49
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggccaagaggcagaattgaagaaacacatctgttaaaagagagacattgtgacaagcaagaactacaaggcagtttgcctggaattgaagccagagccgacaaaacatttgattacaaagcagtttaacaagaaggcggtttaccaagcaggagtgacacaggacc
taagaatgaactcagggaagtgaagagaagagctcaaggagaaaatggaggagataaaacagataaaggatctaattgacaaggattttgataaactcagcaatttggaaattatgaaggaatgcagaaagatatggatgagaagatggacattttaataatacacagaagaactat
aagcttcccttagaagagcaccaaggagcagcagggaactcaggctgatgggaaagactcacagagaaccacagctcaggcccaagaaaatggatggagccagtgaggatcaatggagcaccctgtgtctcttcacaagaagacgatggcaccacaaaaacaaacagggtcact
ggatcccttcacactgtgggacgtgctgcgagaaatgtttgtgtgtctaaagaacaactacaatcgggtatatacgtgcagtgaggactgctgacatctatgtag