

KIF16B-kinesin family member 16B Gene

Reference: PTXBC034984

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	kinesin-like protein KIF16B; C20orf23; KISC20ORF; SNX23; kinesin motor protein; sorting nexin 23; testis secretory sperm-binding protein Li 201a; kinesin family member 16B		

For research use only.

GENE DESCRIPTION	kinesin family member 16B
GENE ID	55614
ORIGIN SPECIES	Human
ACCESSIONS	BC034984
VECTOR	pDONR223
NCBI SYMBOL	KIF16B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgccaggatcaatgcttaccattgaagaagaagtcacaaagacgcttcaggattgcatcgtgtgattagtgaaggctgcagtacatctgcagacacgatgaaggataatgagaaactcacaatggcaccattcaacgtaaactaaatagagcggatggtttctcgtctttggcg
caaatccagatgacctgaaggaccaataaaattagatcccacgctacgtctcgtcgggcaaggaaaggatgcacactcagatttgaggtaagcttgccttgaatttctccaaagaaactatttgaaataaggatgaacgtgtgattgctgagagacgaagtcacttagagaa
atacctcagggaacttttcagcgtgatgctccagtcgcaacatctccctccacatcaacaaagtgggactgactctctgaaacatacattgtgagtttcaccattctcaagaaggagctcttgactacagcagccacgggacgggtag