

SNX10-sorting nexin 10 Gene

Reference: PTXBC034992

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OPTB8; sorting nexin-10; sorting nexin 10		
For research use only.			

GENE DESCRIPTION	sorting nexin 10
GENE ID	29887
ORIGIN SPECIES	Human
ACCESSIONS	BC034992
VECTOR	pDONR223
NCBI SYMBOL	SNX10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttccggaacaacagaaagggaattgttaagtgtctgggttcgagatcctaggatcagaaggaggactctggcattcttactatgagatgtattcatactaatagcatgtgtttacaatgaaacatcctgtgtacgaagaagatagagaattcgtgtggctgaggcaga
gactccaaagtaatgcgttgcgtgtaacaactccagaacttccatctaaaaacctgttttcaacatgaacaatgccagcacgtggatcagcgtcgcagggtctggaagatttctcagaaaagtcctacagaatgcactttgttcagatagcagccttcacctcttctacagagccatct
gaattcagaagacattgagcgtgtgtttctgggcagactaagtactctgtgggaagaagcaattcacaagtttgccttaataatagacgtttccctgaagaagatgaagaaggaaaaaaagaaatgatagattatgattcagaagttcctcctgtggccttgacacagtagtgatgaca
gcagttcacatggatgtaagtaataacagctccgcaggaaatcctga