



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEFF1-transmembrane protein with EGF-like and two follistatin-like domains 1 Gene

Reference: PTXBC035056

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf2; CT120.1; H7365; TR-1; tomoregulin-1; cancer/testis antigen family 120, member 1; transmembrane protein with EGF-like and one follistatin-like domain; transmembrane protein with EGF like and two follistatin like domains 1		

For research use only.

GENE DESCRIPTION	transmembrane protein with EGF-like and two follistatin-like domains 1
GENE ID	8577
ORIGIN SPECIES	Human
ACCESSIONS	BC035056
VECTOR	pDONR223
NCBI SYMBOL	TMEFF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr