

FLJ25770-hypothetical protein FLJ25770 Gene

Reference: PTXBC035224

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ25770
GENE ID	339965
ORIGIN SPECIES	Human
ACCESSIONS	BC035224
VECTOR	pDONR223
NCBI SYMBOL	FLJ25770
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaatcaaaagctgggaatcaataatgaagatctttatcaagtagtggtgtcacatcaatggaggttctcaagttcattttgtgtcatctattcgtgtacaataattgaaaacacatcttcagctgggactttgacacaggttccttttccctaaatgaagtggaaactgatttcctag
aaaaatcatcccatctcctggaagggaacactttgaacgtgttttggagaataatcacatcaagtcgaagatttacagagaagactaaatgaaagcaatgaattgcatgagaacaaaagtgttttattggagcagtcagtcattgatttgcacacaaaacttcaggagatgcaaatggagagag
atgctatggctgacatcagacgaaggagagtcagtcaccaggaggttaagaatcagcttcaaaatacagttcatgaacttgaagctgccaaatgcccttaagaggacatgctgaaagacagcaacacacagatagagcaactacgaaaaatgatccttagtcatgaggagtgcttca
agaaatccggtcaatcctagttgactttgaagaagcctcaggcaaaaaatgtgaacatgacagcatgtctactctgcacttcgcagcttgggtcagctattagtaaaatactaa