

CYorf15B-chromosome Y open reading frame 15B Gene

Reference: PTXBC035312

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	chromosome Y open reading frame 15B
GENE ID	84663
ORIGIN SPECIES	Human
ACCESSIONS	BC035312
VECTOR	pDONR223
NCBI SYMBOL	CYorf15B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaactcacatagcagtagcttcttttcgcgagctttatctgtaggtctgtacttttaacatacatattttggataataaggaggaaacatagcagcagaaaaagaggagggaagaagtcttaagaagtaactgcacatttccaaattactgaaactcaagcccagctggaacagcat
gaaatacacaatgccaaactgcagcaggagaacatggaaatgggagaaaagctaaagagctcactgaccagtagcactgagggaagagtaaacggcgtttgtttttatagacagctgctcaccagtagcaaaagctacacccctaccgtaa