



Reference: PTXBC035359

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD327; CD33L; CD33L1; CD33L2; CDW327; OBBP1; sialic acid-binding Ig-like lectin 6; CD33 antigen-like 1; obesity-binding protein 1; sialic acid binding Ig like lectin 6		

For research use only.

GENE DESCRIPTION	sialic acid binding Ig-like lectin 6
GENE ID	946
ORIGIN SPECIES	Human
ACCESSIONS	BC035359
VECTOR	pDONR223
NCBI SYMBOL	SIGLEC6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]