



Reference: PTXBC035583

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	KIAA0999 protein
GENE ID	23387
ORIGIN SPECIES	Human
ACCESSIONS	BC035583
VECTOR	pDONR223
NCBI SYMBOL	KIAA0999
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]