



Reference: PTXBC035627

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ASC1; HASC-1; asc-1; asc-type amino acid transporter 1; solute carrier family 7 (neutral amino acid transporter light chain, asc system), member 10; solute carrier family 7, (cationic amino acid transporter, y+ system) member 10; solute carrier family 7.		

For research use only.

GENE DESCRIPTION	solute carrier family 7, (neutral amino acid transporter, y+ system) member 10
GENE ID	56301
ORIGIN SPECIES	Human
ACCESSIONS	BC035627
VECTOR	pDONR223
NCBI SYMBOL	SLC7A10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]