



ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

Reference: PTXBC035683

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	serine/threonine protein kinase PIK3CG; PI3CG; PI3K; PI3Kgamma; PIK3; p110gamma; p120-PI3K; phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit gamma isoform; 1-phosphatidylinositol 3-kinase; PI3-kinase subunit gamma; phosphatidylinositol 3 k		

GENE DESCRIPTION	phosphoinositide-3-kinase, catalytic, gamma polypeptide
GENE ID	5294
ORIGIN SPECIES	Human
ACCESSIONS	BC035683
VECTOR	pDONR223
NCBI SYMBOL	PIK3CG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1