



Reference: PTXBC035718

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EBAF; LEFTA; LEFTYA; TGFB4; left-right determination factor 2; TGF-beta-4; endometrial bleeding-associated factor; left-right determination factor A; protein lefty-2; protein lefty-A; transforming growth factor beta-4		

For research use only.

GENE DESCRIPTION	left-right determination factor 2
GENE ID	7044
ORIGIN SPECIES	Human
ACCESSIONS	BC035718
VECTOR	pDONR223
NCBI SYMBOL	LEFTY2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]