



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C12orf53-chromosome 12 open reading frame 53 Gene

Reference: PTXBC035736

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf53; LEDA1; PANP; leda-1; PILR alpha-associated neural protein; PILR-associating neural protein; liver endothelial differentiation-associated protein-1; paired immunoglobulin-like type 2 receptor-associating neural protein; PILR alpha associated neural		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 53
GENE ID	196500
ORIGIN SPECIES	Human
ACCESSIONS	BC035736
VECTOR	pDONR223
NCBI SYMBOL	C12orf53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagtcaggagatgtgctgctgcgtctgctgtcccaactctctctctctgccaactgctgttctgctgccctccaccgcctgctcagggctctctcatctccctcgaacccaccagccccagcccccctgtgcaggggagagccctcgccccacgtcatgtgtgcgtgtgggagcgcagcactcccaacgcgcatctctcgggtcccaagatcacgtcgtgccaagtctgctgctgcactgcacccccagccacccccatcaggctttgaggagggggccgcctcatccaataacccctgggctatcgtgtgggtgccaccgtgtctcgcagaggaatgagggggacc
ccaactctgccaatcccggaattctgcaatctgttttgcagccctcatgggtcgcacaacccaaccccaactcagactcgaatcgagagtgatggagatggcttactgtgagagagcactgccacctgcagccatctgttcggggccgtggggaaggtgtggacccccact
ctatgtcaaatattacatctcatcatcatctgttcggcactgcacatctcaagtctcgtggagcagcaggaagcgcgcagcagccctcaggcgcaggtgtgagcgttcgtgagcagggagcagcagccactgacagactgtccccggctgagtcactgtgtcgtgg
ggcctctggagctaccactccccccccctgacatcagcagggcccccagggggcccccgcggccctggatggcccccaacccaagggtccagcttcagatgaagcattccccctgtggaatctgtga

More info: Please contact us for any question: contact@proteogenix.fr