



Reference: PTXBC035750

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CYSLT1; CYSLT1R; CYSLTR; HMTMF81; cysteinyl leukotriene receptor 1; G-protein coupled receptor HG55; LTD4 receptor; cysteinyl leukotriene D4 receptor		

For research use only.

GENE DESCRIPTION	cysteinyl leukotriene receptor 1
GENE ID	10800
ORIGIN SPECIES	Human
ACCESSIONS	BC035750
VECTOR	pDONR223
NCBI SYMBOL	CYSLTR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1