



Reference: PTXBC035763

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EXO84; Exo84p; SEC84; exocyst complex component 8; exocyst complex 84 kDa subunit		

For research use only.

GENE DESCRIPTION	exocyst complex component 8
GENE ID	149371
ORIGIN SPECIES	Human
ACCESSIONS	BC035763
VECTOR	pDONR223
NCBI SYMBOL	EXOC8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr