



Reference: PTXBC036029

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	metalloproteinase-disintegrin ADAM22-3; ADAM 22; MDC2; disintegrin and metalloproteinase domain-containing protein 22; a disintegrin and metalloproteinase domain 22; metalloproteinase-like, disintegrin-like, and cysteine-rich protein 2; ADAM metallopeptid		

GENE DESCRIPTION	ADAM metallopeptidase domain 22
GENE ID	53616
ORIGIN SPECIES	Human
ACCESSIONS	BC036029
VECTOR	pDONR223
NCBI SYMBOL	ADAM22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]