



Reference: PTXBC036077

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERF3B; GST2; eukaryotic peptide chain release factor GTP-binding subunit ERF3B; eukaryotic peptide chain release factor subunit 3b; G1 to S phase transition 2		

For research use only.

GENE DESCRIPTION	G1 to S phase transition 2
GENE ID	23708
ORIGIN SPECIES	Human
ACCESSIONS	BC036077
VECTOR	pDONR223
NCBI SYMBOL	GSPT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]