



Reference: PTXBC036101

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Fuc-TIX; alpha-(1,3)-fucosyltransferase 9; fucT-IX; fucosyltransferase 9 (alpha (1,3) fucosyltransferase); fucosyltransferase IX; galactoside 3-L-fucosyltransferase; fucosyltransferase 9		

For research use only.

GENE DESCRIPTION	fucosyltransferase 9 (alpha (1,3) fucosyltransferase)
GENE ID	10690
ORIGIN SPECIES	Human
ACCESSIONS	BC036101
VECTOR	pDONR223
NCBI SYMBOL	FUT9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]