



Reference: PTXBC036113

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM6A; PRLMNS; hDIS3L2; DIS3-like exonuclease 2; DIS3 mitotic control homolog-like 2; family with sequence similarity 6, member A; DIS3 like 3'-5' exoribonuclease 2		

For research use only.

GENE DESCRIPTION	DIS3 mitotic control homolog (S. cerevisiae)-like 2
GENE ID	129563
ORIGIN SPECIES	Human
ACCESSIONS	BC036113
VECTOR	pDONR223
NCBI SYMBOL	DIS3L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]