

C1orf71-chromosome 1 open reading frame 71 Gene

Reference: PTXBC036200

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf71; PPP1R64; consorin; protein phosphatase 1, regulatory subunit 64; consorin, connexin sorting protein		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 71
GENE ID	163882
ORIGIN SPECIES	Human
ACCESSIONS	BC036200
VECTOR	pDONR223
NCBI SYMBOL	C1orf71
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgacagcgatactcctacatattatctgcaaatagaaccacaagatggatgtcatctgtgacagcgtggaaaggagtgacctgtctgcttctgcatcagatgaaatgaaatcagcttgacggggacgggcatgagcatctgaccagcagtgacagtgcatgggaaagc
cccaagtgtctgagcaggacagtcataataatgaaagctgcacattgagctgcgaggtggctgagtgagaaacttgcaaaacacccctttgtgaagcctccagagatgaacaggcctcttgggaaaggacaaaaaattcctggaaaaagaagccaagaagcaaaaaaggact
gctaagaagataccaccaggactttttcaggagatattgcacctttaatgcaagaaaaagtactaagcgcagtcacatatgctgttgatgatgaagaagctgctgaagttaattgtaagagcagccagaggcgccaaagctgttctgagctctgttttactatagaggtgaagtgag
cagttggattcaagagcacttccctttgccttcacatagcagaatcctatttcaggaggagcagatgagaaagcaatgaaattcattcagctagaacgattgtatcatgagcaattgctcgcaaatcttctgccattcaagaacagtgagggttaagtacagccaacatggaatgtgaa
ttagcagtaatccaacaagtcactatctacaaaattaacatattgacctggtttataa