



Reference: PTXBC036216

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CUL-4B; MRXHF2; MRXS15; MRXSC; SFM2; cullin-4B; cullin 4B		
For research use only.			

GENE DESCRIPTION	cullin 4B
GENE ID	8450
ORIGIN SPECIES	Human
ACCESSIONS	BC036216
VECTOR	pDONR223
NCBI SYMBOL	CUL4B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr