

C20orf12-chromosome 20 open reading frame 12 Gene

Reference: PTXBC036231

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat-containing protein C20orf12; C20orf12; ANKRD64; C20orf84; bA189K21.1; bA189K21.8; dJ568F9.2; double zinc ribbon and ankyrin repeat-containing protein 1; ankyrin repeat domain 64; double zinc ribbon and ankyrin repeat domains 1		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 12
GENE ID	55184
ORIGIN SPECIES	Human
ACCESSIONS	BC036231
VECTOR	pDONR223
NCBI SYMBOL	C20orf12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgctggttcagtggtgtccctcagatcaccactacagtgctcagcctggaagctaacatgaaatgataacaatcgttttgaaatgaatcagacactccagatgtcaacatattatactctggatggcagcaaacctgaatttctaaagagaattggtatggggaaa
ataacacatttaagtataataaacctattactctgcctgatggaaaaatacaagttaaagctattgtgtctctaaagactgcagacagagtggtgcatggcacaaggtgtttcacgttagactatgaaccaccaatagctctcctgaagacaattgtgaaatgttctcaagattctcaaggc
aggaattcaaaatgattgttgatcaaaactaaagaaaaataaagaactctgaaatcaacgcagctgggaatgttaaccttagaaagtctccagtcctcccggttttgcacacgttaagcggtcagaagtgttgacaagcacggagataatgagaattcaagggaagacagacttctca
agtgtgccactgcctgcccccgccatcagatcccttgcctcgtctgtcaagaatgtggtctctgtccacccatattggtgtgtctccacccccagaaggagctcagatgggtgtgtgtcagaatgcagaagcttggtacctgaacactccatctgcgtgtgtgtga