

RP4-691N24.1-ninein-like Gene

Reference: PTXBC036380

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NLP; ninein-like protein; ninein like		
For research use only.			

GENE DESCRIPTION	ninein-like
GENE ID	22981
ORIGIN SPECIES	Human
ACCESSIONS	BC036380
VECTOR	pDONR223
NCBI SYMBOL	RP4-691N24.1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr