



Reference: PTXBC036444

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cytoplasmic polyadenylation element-binding protein 3; CPE-BP3; CPE-binding protein 3; hCPEB-3; cytoplasmic polyadenylation element binding protein 3		

For research use only.

GENE DESCRIPTION	cytoplasmic polyadenylation element binding protein 3
GENE ID	22849
ORIGIN SPECIES	Human
ACCESSIONS	BC036444
VECTOR	pDONR223
NCBI SYMBOL	CPEB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr