

RYBP-RING1 and YY1 binding protein Gene

Reference: PTXBC036459

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ring1 interactor RYBP; AAPI; APAP-1; DEDAF; YEAF1; RING1 and YY1-binding protein; DED-associated factor; YY1 and E4TF1 associated factor 1; apoptin-associated protein 1; death effector domain-associated factor; RING1 and YY1 binding protein		

For research use only.

GENE DESCRIPTION	RING1 and YY1 binding protein
GENE ID	23429
ORIGIN SPECIES	Human
ACCESSIONS	BC036459
VECTOR	pDONR223
NCBI SYMBOL	RYBP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccatggcgacaagaagagcccgaccaggcccaaaagacaagcgaacacctgccacagacgaagggtttgggattgtagcgtctgcaccttcagaacacgtgctgaagcctttaatgcagcatctgcgatgtgaggaaaggcacctccaccagaaaacctcggatcaattctca
gctgggtggcacaacaagtggcacaacagatagccacccaccacccctctaaaaggagaagaaggagaaggttgaagcaggacaaagagaacacgtgagaaagacaaggaaattagtcctagtgtttaccaagaaaaataccaacaagaaaacaaacaaagtctgacattctgaa
agatcctcctagtgaagcaaacagcacaagctctgcaaatgctacaacaagaccagcgaacaaatcacacctcaaggcccccgtgaaaaacgtggacaggagcactgcacagcagttggcagtaactgtgggcaacgtcacctgcatcagacactttaagaaaagactcgt
cctcatcgacatcctcatccacagtgacctccagtgcagggtcagaacagcagaaccagagcagctcgggtcagagagcacagacaagggtcctccctcctccacgccaaggcgacatgtcagcagtcgaatgatatttttga