



Reference: PTXBC036465

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRESEK; IFAP; KFSD; KFSDX; OLMSX; S2P; membrane-bound transcription factor site-2 protease; SREBPs intramembrane protease; endopeptidase S2P; keratosis follicularis spinulosa decalvans; membrane-bound transcription factor protease, site 2; site-2 protease		

For research use only.

GENE DESCRIPTION	membrane-bound transcription factor peptidase, site 2
GENE ID	51360
ORIGIN SPECIES	Human
ACCESSIONS	BC036465
VECTOR	pDONR223
NCBI SYMBOL	MBTPS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]