



Reference: PTXBC036515

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3-C5; METH1; A disintegrin and metalloproteinase with thrombospondin motifs 1; ADAM-TS 1; ADAM-TS1; ADAMTS-1; METH-1; a disintegrin-like and metalloprotease (reprolysin type) with thrombospondin type 1 motif, 1; human metalloproteinase with thrombospondi		

GENE DESCRIPTION	ADAM metallopeptidase with thrombospondin type 1 motif, 1
GENE ID	9510
ORIGIN SPECIES	Human
ACCESSIONS	BC036515
VECTOR	pDONR223
NCBI SYMBOL	ADAMTS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]