

COL1A1-collagen, type I, alpha 1 Gene

Reference: PTXBC036531

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EDSC; OI1; OI2; OI3; OI4; collagen alpha-1(I) chain; alpha-1 type I collagen; alpha1(I) procollagen; collagen alpha 1 chain type I; collagen alpha-1(I) chain preproprotein; collagen of skin, tendon and bone, alpha-1 chain; collagen, type I, alpha 1; pro-a		

For research use only.

GENE DESCRIPTION	collagen, type I, alpha 1
GENE ID	1277
ORIGIN SPECIES	Human
ACCESSIONS	BC036531
VECTOR	pDONR223
NCBI SYMBOL	COL1A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr