

C13orf36-chromosome 13 open reading frame 36 Gene

Reference: PTXBC036540

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf36; serine-rich and transmembrane domain-containing protein 1; serine rich and transmembrane domain containing 1		

For research use only.

GENE DESCRIPTION	chromosome 13 open reading frame 36
GENE ID	400120
ORIGIN SPECIES	Human
ACCESSIONS	BC036540
VECTOR	pDONR223
NCBI SYMBOL	C13orf36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgaacctgacacttctcaggatttctgggaagtgtggagaatggaactttcttgagctgttccacatccctgtccacgtcagtgaccatctcagccacctgtcaaacgtctacatctatgtgccatattcctcagccttttagcgtttctgtctctgttttaacattgcctccag
aggctcaaaaatatcatctctcagttctctaccagggtatccaagcgacgtggaagtctttcaccaatttgaagctcgcagcatttctctcagaggtccacttttcaaacctttcatctga