



Reference: PTXBC036661

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXCR7; CMKOR1; CXC-R7; CXCR-7; GPR159; RDC-1; atypical chemokine receptor 3; C-X-C chemokine receptor type 7; G protein-coupled receptor; G-protein coupled receptor 159; G-protein coupled receptor RDC1 homolog; chemokine (C-X-C motif) receptor 7; chemokine		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) receptor 7
GENE ID	57007
ORIGIN SPECIES	Human
ACCESSIONS	BC036661
VECTOR	pDONR223
NCBI SYMBOL	CXCR7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]