

LRRC18-leucine rich repeat containing 18 Gene

Reference: PTXBC036722

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UNQ933; UNQ9338; VKGE9338; leucine-rich repeat-containing protein 18; leucine rich repeat containing 18		

For research use only.

GENE DESCRIPTION	leucine rich repeat containing 18
GENE ID	474354
ORIGIN SPECIES	Human
ACCESSIONS	BC036722
VECTOR	pDONR223
NCBI SYMBOL	LRRC18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttaagggtgagaaagcccaaggcgaagaagatcaccctcaagggtggccagggaattgcatcaaaatcattttgatgggaaaaagcgcttgacttgagcaagatgggaattaccacctcccaagtgtattctgcgccttagtgacatggacgagctggaccttagccggaatct
tatcaggaagatccctgactccatctccaagttccagaacctccgggtgctggacctgcacagcaactacatagacaagctgcctgagtcattggccagatgaccagcctgtctacctcaacgtcagcaacaacccgctgaccagcaacggctgcccgtggagctgaagcaactcaa
gaacatccgctgtgaacctaggctgaaccacctggacagcgtgtccaccacactggggccctgaaggagctccacaggttagggctcatgacaacctactgaacaacatccccgtgagcatctccaagctcccaagctgaaaaagctcaacataaagggaaccccttccaa
agccaggtgagtcggaatatcatagactccatcaggaggctggagaactgtgtgttggaggagaaggatctgtgtcgccttgagaaaatgccaaaacgcccgggacaacctgaatagaatcaagaacatggccacgacgacaccgagaaaagaccatcttccaatctga
tctcaccatccatgccaaggactcctgggaagactggagaatacgttgacatcttccag