

FCGR3A-Fc fragment of IgG, low affinity IIIa, receptor (CD16a) Gene

Reference: PTXBC036723

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD16A; FCG3; FCGR3; FCGR3III; FCR-10; FCRIII; FCRIIIA; IGFR3; IMD20; low affinity immunoglobulin gamma Fc region receptor III-A; CD16a antigen; Fc fragment of IgG low affinity IIIa receptor; Fc fragment of IgG, low affinity III, receptor for (CD16); Fc fra		

For research use only.

GENE DESCRIPTION	Fc fragment of IgG, low affinity IIIa, receptor (CD16a)
GENE ID	2214
ORIGIN SPECIES	Human
ACCESSIONS	BC036723
VECTOR	pDONR223
NCBI SYMBOL	FCGR3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtggaggggctgggaaaggctgtttacttctctgtctagtgcgtttgtcccttttagggctccggatattctgtgactgtctctccagtggtggcatatgtggcagctgctctcccaactgctctgctactttagttcagctggcatgaggactgaagatctccaaaggctgt
gggttctctggagcctcaatggtacaggggtgctcgagaaggacagtggtgacttgaagtgcaggagcctactccccctgaggacaattccacacagtggtttcacaatgagagcctcatctcaagccaggcctcgagctactcattgacgctgccacagtcgacgacagtgagagtag
agggtccagacaaacctctccaccctcagtgacccgggtgcagctagaagtcacatcggtgctgtgtctccaggccccctcgggtgggtttcaaggagggaagacccattacacgtgaggtgtcacagctgggaagaacactgctctgataaggtcacatattacagaatggcgaaggca
ggagatttttcatcataattctgacttcatattcaaaagccacactcaagacagcggctcctacttgcagggggctgtgtggagtaaaaatgtgtcttcagagactgtgaacatcaccatcactcaagggtttggcagtgtaacctatcatcattcttcacctgggtaccagctcttt
ctgcttgggtgatgtactccttttgcagtggtgacacaggactatattctctgtgagacaacattcgaagctcaacaagagactggaaggaccataaatttaattggaagaaggacccctcaagacaaatga